

SuPrime HF Premix (2X)

Product Name	Cat. No.	Size
SuPrime HF Premix (2X)	HF-2000	1.0 ml X 1


Package information

HF-2000	2X SuPrime HF Premix (1.0 ml X 1) - with SuPrime HF DNA Polymerase, dNTPs mixture, reaction buffer, enzyme stabilizer and loading dye
---------	--

Description

SuPrime HF Premix contains SuPrime HF DNA Polymerase, dNTPs mixture, reaction buffer, enzyme stabilizer and loading dye for electrophoresis.

〈Error Rate of SuPrime HF DNA Polymerase〉



► Genomic DNA were amplified with SuPrime HF DNA Polymerase and other Polymerase. 1 kb PCR products were cloned into vector. Each 100 clones were selected and subjected to sequence analysis to determine the error rate.

Usage Information

- SuPrime HF Premix produce **blunt end DNA fragments**.
- The extension time for long PCR is **20~30 sec/kb**.
- The denaturation and extension temp. is **98℃** and **68℃**.
- The concentration of SuPrime HF Premix is **2X**.
- The amount of SuPrime HF DNA Polymerase per 10 μ l SuPrime HF Premix is **1 unit**.
- If the smearing or non-specific products are appeared, use the SuPrime HF DNA Polymerase (Cat. No. HF-1000).

● **Research Use Only**

● **Store at -20℃**

Protocol

The following 20 μ l reaction volume can be used for PCR.

1. Prepare the following components to a PCR tube.

Components	Volume
DW	add up to 20 μ l
2X SuPrime HF Premix	10 μ l
Upstream Primer (10 pmoles/ μ l)	0.5 ~ 1 μ l
Downstream Primer (10 pmoles/ μ l)	0.5 ~ 1 μ l
Template DNA *	X μ l

* Amount of template DNA

- Plasmid, Lambda DNA, BAC DNA: 1 pg~5 ng
- Genomic DNA: 10 ng~250 ng

2. PCR cycling

Step	2-step PCR		3-step PCR		Cycles
	Temp.	Time	Temp.	Time	
Initial denaturation	98℃	2 min	98℃	2 min	1
Denaturation	98℃	10 sec	98℃	10 sec	25~35
Annealing	-	-	X℃	20 sec	
Extension	68℃	10~30s/kb	68℃	10~30s/kb	
Final Extension	68℃	5 min	68℃	5 min	1

3. Separate the PCR products by agarose gel electrophoresis and visualize with EtBr or any other means.