

PrimePrep™ Gel Purification Kit

[지금 받으신 제품은 Sample size로써, 10회 prep.을 하실 수 있습니다.
Buffer PW에 첨가할 Ethanol의 양은 bottle label에 표기되어있습니다.]

Introduction

PrimePrep™ Gel Purification Kit offer simple, rapid and cost-effective method for purification from agarose gel in TAE or TBE buffer system..
The purified DNA can be directly used in ligation, sequencing and other downstream application.

Kit Components

Cat. No. Reagents	K-8000 (50 prep.)
Spin column	50 ea
Buffer Gel-B	60 ml (30 ml x 2)
Buffer PW	10 ml
Buffer PE	10 ml

Before you begin

- Add ethanol to Buffer PW before use.
→ bottle의 label에 첨가량이 표시되어 있습니다.

Experimental Protocol

1. Cut out the interesting DNA fragment with a sharp scalpel or razor blade.
Minimize gel volume by cutting gel slice as small as possible.
2. Weigh the gel slice in a microcentrifuge tube. Add 3 volumes of Buffer Gel-B to 1 volume of gel.
If the 100 mg of agarose gel slice, add 300 ul of Buffer Gel-B.
If more than 1.5% agarose gel, add 6 volumes of Buffer Gel-B.
3. Vortex the mixture and incubate at 55°C until the agarose gel is completely melted (5~10min). To help the dissolving of gel, vortex the tube every 2~3 min during the incubation. Centrifuge the tube briefly at room temperature.
4. (Optional) Add 1 gel volume of isopropanol to the sample and vortex to mix.
Do not centrifuge after mixing well.
This step increases the yield of DNA fragment <200bp or >5kb.
5. Transfer the mixture to a Spin column.
6. Centrifuge for 1 min at 13,000 rpm. Discard the flow-through and re-inserting the spin column to the collection tube.
7. (Optional) Add 500 ul of Buffer Gel-B to the column and centrifuge for 30 sec at 13,000 rpm. Discard the flow-through and re-inserting the spin column to the collection tube.
This step is complete removal of agarose for very sensitive applications.
8. Add 700 ul Buffer PW and centrifuge for 30 sec. at 13,000 rpm.
Discard the flow-through and re-inserting the spin column to the collection tube.
9. Centrifuge for an additional 1 min at 13,000 rpm to remove residual wash buffer.
Residual ethanol of washing buffer may inhibitor subsequent enzymatic reaction.
10. Transfer the spin column to new 1.5 ml microcentrifuge tube.
The 1.5 ml microcentrifuge tube is not provided.
11. Add 50 ul of Buffer PE or deionized distilled water to the center of the membrane in the column, let stand for 1 min and centrifuge for 1 min at 13,000 rpm.
For larger fragment(>5kb), use pre-warmed (70°C) Buffer PE for best efficiency.