PrimePrep[™] Viral RNA/DNA Extraction Kit

Introduction

PrimePrepTM Viral RNA/DNA Extraction Kit is provide a fast, easy method for the preparation of viral RNA and DNA from plasma, serum, cell-free bo dy fluids, cell-culture supernatants and virus-infected samples. *PrimePrep*TM Viral RNA/DNA Extraction Kit buffer system provides the effective binding c ondition of RNA and DNA to microfiber-silica-based membrane through m ix with lysis and binding buffers. And then the impurities on the membrane are washed away by two different wash buffers. Purified viral RNA/DNA re ady for use in downstream applications such as PCR, RT-PCR, cDNA synthe sis and real-time PCR etc.

Kit Components

Cat. No. Reagents	KR-2000 (50 Prep.)
Spin column	50 ea
Collection tube	100 ea
Buffer VRL	40 ml
Buffer VRW1	20 ml
Buffer VRW2	11 ml
Buffer VRE	10 ml

Before you begin

- Add 20 ml absolute ethanol to Buffer VRW1 before use.
- Add 44 ml absolute ethanol to Buffer VRW2 before use.
- Prepare β-mercaptoethanol (14.2 M), isopropanol, 1.5 ml microcentrifuge tube.

! Caution:

- This product is for research use only.
- Buffer VRL and Buffer VRW1 contains strong denaturant. Be careful to avoid contacting with skin and eyes. In the case of such contact, wash immediately with plenty water.



Experimental Protocol

- 1. Prepare 200 µl of sample (plasma, serum, cell-free body fluids, cell-culture supernatants, virus-infected samples) into a 1.5 ml microcentrifuge tube.
- If the sample volume is larger than 200 μ l, increase the amount of Buffer VRL (e.g. a 400 μ l sample will require 700 μ l Buffer VRL).
- **VNote**: Virus-infected samples are infected feces and cell line. In the case of feces samples, transfer 20~50 mg in 1.5 ml microcentrifuge tube and dissolve with 250 μ l of DEPC-DW and vortex for 10~15 sec. Centrifuge at 13,000 rpm for 30 sec at room temperature. Use approximately 150~200 μ l of supernatant.
- 2. Add 350 μ l of Buffer VRL and then add 3.5 μ l β -mercaptoethanol (β -ME, 14.2 M). Mix by vortexing for 10~15 sec.
 - \cdot Check Buffer VRL for precipitation. If happened, dissolve precipitate completely by incubation at 40~50 °C until dissolve.
 - \cdot $\beta\text{-ME}$ must be add to Buffer VRL as ratio of 10 μI of $\beta\text{-ME}$ per 1 ml of Buffer VRL.
- 3. Incubate at room temperature (15~25 °C) for 10 min. After this step, cen trifuge at 3,000 rpm for 5 sec at room temperature.
 - For protect from contamination by carryover, spin-down the tube to remove dr ops from the inside of the lid.
- 4. Add 150 μ of isopropanol to lysate, and mix by pulse-vortexing for 15 sec. After this step, briefly centrifuge the tube to remove drops from the inside of the lid. Transfer the mixture to spin column (fit in a 2 ml collection tube).
- 5. Centrifuge at 13,000 rpm for 30 sec at room temperature. And discard the pass-through and reinsert the spin column into the same collection tube.
- 6. Add 700 μ l of Buffer VRW1 to the spin column. Centrifuge at 13,000 rpm f or 30 sec at room temperature. And discard the pass-through and reinsert t he spin column into the new collection tube.
- 7. Add 700 μ l of Buffer VRW2 to the Spin column. Centrifuge at 13,000 rpm f or 30 sec at room temperature. And discard the pass-through and reinsert t he spin column into the new collection tube.
- 8. Centrifuge at 13,000 rpm for 1 min at room temperature to remove residual wash buffer. Transfer the spin column to a new collection tube.
- Residual ethanol may interfere with downstream reactions.
- 9. Add 50~150 μ l of Buffer VRE to the center of the membrane in spin Col umn, and wait for at least 1 min at room temperature.
- 10. Centrifuge at 13,000 rpm for 1 min at room temperature.
 - \cdot Purified RNA/DNA can be stored at -20 $^\circ\!C$ for immediate use and stored at -70 $^\circ\!C$ for long term storage.

