

## 100 bp Plus DNA Marker

| Product Name           | Cat. No. | Size       |
|------------------------|----------|------------|
| 100 bp Plus DNA Marker | M-2000   | 0.5 ml X 1 |
|                        | M-2002   | 0.5 ml X 5 |

### Specification


Concentration: 83 ng/ $\mu$ l

Recommended loading: 4.0 ~ 5.0  $\mu$ l/5 mm lane width

Number of band: 18

Size range: 100 ~ 10,000 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,  
 5% Glycerol, 0.005% Bromophenol Blue,  
 0.005% Xylene Cyanol



5.0  $\mu$ l loading/ 5 mm lane width  
 1.0% TAE agarose gel with stained EtBr

- Research Use Only
- Store at -20°C