

100 bp DNA Marker

Product Name	Cat. No.	Size
100 bp DNA Marker	M-1000	0.5 ml X 1

Specification

Concentration: 66 ng/ μ l

Recommended loading: 3.0 μ l/ 5 mm lane width


Number of band: 11

Size range: 100 - 1,500 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,
5% Glycerol, 0.005% Bromophenol Blue,
0.005% Xylene Cyanol

☐ Research Use Only

☐ Store at -20°C



3.0 μ l/ 5 mm lane width

1.5% TBE agarose gel with stained EtBr