

100 bp Plus DNA Marker

Product Name	Cat. No.	Size
100 bp Plus DNA Marker	M-2000	0.5 ml X 1

Specification


Concentration: 80.7 ng/ μ l

Recommended loading: 3.0 ~ 4.0 μ l/5 mm lane width

Number of band: 19

Size range: 100 ~ 10,200 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,
2.5% Ficoll, 0.005% Bromophenol Blue,
0.005% Xylene Cyanol



● **Research Use Only**

● **Store at -20°C**

3.0 μ l loading/ 5 mm lane width
1.0% TBE agarose gel with stained EtBr